

# Guiding Principles

- 1. The reagents core should serve the entire community.  
(It is possible to serve the community and to advance diagnostics.)**
- 2. Validation, annotation, coordination, and curation into a database are the missing links, and the most valuable service the Reagents Core can provide.**
- 3. New, higher-throughput technologies should be developed in preparation for expanding this project to target the entire proteome.**

## Sources

## Applications

## Validation/Annotation

commercially available

unscreened hybridomas

requests from community

antibodies funded by grants




**In what types of applications should  
affinity reagents be required to  
perform?**

# Research Areas of Antibody Customers

- 56% on cellular development
- 47% on gene/protein regulation
- 44% on molecular mechanisms of human diseases
- 39% on molecular structure/function characterization
- 28% on non-human disease

## CST Market analysis



### Shared by Abnova

- WB
- IHC
- IP
- ELISA

### Shared by Epitomics

- 79% Western
- 60% ELISA
- 55% Immunofluorescence
- 50% IP
- 44% IHC

### Shared by Neoclone

Survey of most recent 100 projects


- 58% Western
- 34% ELISA
- 28% IHC
- 9% Immunoaffinity Purification
- 6% IP
- 6% FACS
- 2% Chromatin IP (ChIP)

### Antibodies by Design, a Division of MorphoSys

- Westerns
- IHC
- ELISA
- FACS
- IP

# SISCAPA: an emerging technology.

Journal of Proteome Research 2004, 3, 235-244



>10<sup>4</sup> fold enrichment S2N

# Recommended Target Applications

## Sources

commercially available

requests from community

unscreened hybridomas

## Applications

1. Western
2. ELISA
3. Immunofluorescence / FACS
4. Immunoprecipitation
5. Immunohistochemistry
6. SISCAPA

## Validation/Annotation

application-independent

application-dependent

# Priority Questions

1. How sensitive is the antibody (affinity/avidity)?

- Biacore?
- Antigen spike experiments?

2. Does my antibody bind what I want it to bind?

- tissue arrays (+/- controls)
- protein or prEST arrays
- cell lines +/- expression of target
  - siRNA
  - cDNA
  - KO

3. How specific is my antibody? (S2N)

- tissue arrays (+/- controls)
- protein or prEST arrays
- cell lines (+/- controls)
- X-SPECIES??



# Validation methods for antibodies (application-independent)

Method	Description	Examples
Antigen-based	Assays based on the antigen used for immunization (immunogen)	ELISA, protein arrays, SPR, antigen adsorption
Target-based	Analysis of native or partially denaturated protein from natural sources (such as cell lysates)	Western blots, IHC, immunocapture (pull-downs)
RNA-based	Comparison of expression levels at the protein and RNA level	Transcript profiling, in situ hybridizations
DNA-based	Bioinformatic analysis using predictive algorithms (as compared to experimental data)	Signal peptide, transmembrane regions, localization signals
Genetics-based	The use of genetic mutants or recombinant constructions to validate the target	Transgenetics, RNAi, GFP-fusions (subcellular localization)
Epitope-based	Comparison of two or more antibodies directed to different parts of the same target	Antibodies to PrESTs or synthetic peptides

Uhlen et al (2005) A human protein atlas for normal and cancer tissues, Mol Cell Proteomics, in press (on-line)

# **Application-independent**

- 1. What are the binding affinities of the Ab? Include Biacore measurements with each Ab.**
- 2. Characterize specificity by screening on arrays of expressed human proteins (prESTs?), reporting prevalence of cross-reactions.**
- 3. Validate target identity with control lysates (siRNA?).**
- 4. Better information on species cross-reactivity (human, mouse, rat, etc)**
- 5. Include specific antigen information/sequence**
- 6. Formulation of Abs (carriers and buffers included in mix)**
- 7. Accurate shelf life and storage conditions**
- 8. Provide both Ab concentration (mg/ml) and working dilutions (1:1000 etc.)**
- 9. Showing enough data from reagents is always an issue. If groups could share more data with the companies without publication issues, then perhaps it could benefit everyone using the reagents.**
- 10. Test the new lots in all the applications side by side with the old lot.**

# Application-Specific Validation

1. **Western blot**: show entire blot with positive and negative controls, show background signal, results on multiple cell lines, cross-species “zoo blots”
2. **ELISA**: standard curve showing full dynamic range (upper and lower limits of detection) and c.v. of analyte spiked into plasma/cell lysate.
3. **Immunofluorescence**: confocal images from multiple cell lines and tissues, including positive and negative controls where possible.
4. **Immunoprecipitation**: western blots of both initial lysate and captured proteins, show WB data for protein complex pulldowns
5. **SISCAPA**: MS-based determination of enrichment, specificity, and limit of quantification of target peptide relative to isotopically labeled peptide spiked into complex mixture (plasma or cell lysate).

# Immunohistochemistry

Formalin-fixed, paraffin-embedded tissue specimens afford the best morphology and represent the broadest possible range of tissues and diseases.

- Screen tissue arrays and cell line array (positive and negative control cell lines in paraffin blocks)
- For ligands that have been studied previously, the pattern and anatomic distribution of staining should be consistent with what has been reported in the literature.  
(immunohistochemistry, Northern-blot analysis, in situ hybridization, EST expression profiling, radioimmunoassay, and ligand-binding data)
- For “orphan ligands,” supporting data are not available, and evaluation of the specificity and validity of IHC findings relies more heavily on cross-concordant findings with multiple antibodies and on gene expression studies such as Northern-blot analysis or in situ hybridization.
- Another criterion for evaluating antibody specificity and performance is that the subcellular localization of the staining signal should reflect the protein's function.